

Alcohol, ketone and aldehyde measurements during TRACE-P

Eric Apel, Daniel Riemer

Richard Lueb

Alan Hills

Sebastian Zindel

Scott Eisele

Instrument: GC/MS
Time resolution: 5 minutes
Sensitivity: 5-30 pptv

Compounds measured:

Alcohols

MeOH

EtOH

Ketones

Acetone

Butanone

Pentanone

Aldehydes

Acetaldehyde

Propanal

Butanal

Pentanal

DATA STATUS

As of latest update (Nov. 9) -

Alcohols

Methanol and Ethanol - OK


Ketones

Acetone and Methyl Ethyl Ketone - OK

2-pentanone - More work needed

Aldehydes

Acetaldehyde, Propanal, Butanal and Pentanal - More work needed



RACK #51

FOCUS
RACK DIVISION
SPHERIC RESEARCH
RACK

RACK #51
EMPTY

SCOPTRE



Output file:

Add text here, enter with F8

operation in laboratory

bypass A/D cards	bypass Wallow	bypass ICATS	GC 2	close aircraft inlet	write data file
<input type="checkbox"/>					
Sample load	80	GC 1	auto	GC1	
Cryofocus trap cool	30	T1	5	GC2	
pre cryo trap flush	30	Ramp	70		
Cryofocus	90	T2	5		
Injection	5	Cooldown	10		

NCAR--Fast GC/MS System



static air T (C)	0.0	solar elev. relative to DC-8	0.0
aerosols (particles/ml)	0	CH4 (ppm)	0.00
GPS altitude	0	cryogenic FP (C)	0.0

Mode: **In-flight run** Exit

GMT time: latitude:
 longitude:



load sample	cool cryo trap	pre-cryo trap flush	cryofocus	inject	T1	T ramp	T2	cooldown GC
<input type="checkbox"/>								

Gas Flow

sample volume (ml)	87.5	He carrier (sccm)	3.85	grab sampling	<input type="checkbox"/>
sample (sccm)	-0.9	He backflush (sccm)	0.06	grab sample flow (slpm)	2.53
water trap cooling (slpm)	-0.02	water trap flush (sccm)	0.0	grab sample duration (s)	36
sample trap cooling (slpm)	0.22	BPR (Torr)	801		
cryo trap cooling (slpm)	0.02		800		


GC/MS runs:






GC 1 Fan/shutter	MS start	GC 2 Fan/shutter
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
sample trap flush	water trap flush	column flush
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Valco 1	Valco 2	Valco 3
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Sample	Aircraft inlet	Pump isolation
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Temperature (C)

Probe tip	1371.0	Water trap	35.3
	35.0		-10.0
Probe inlet	1371.0	Sample trap heat	74.1
	35.0		-110.0
Inlet line	1371.0	Sample trap cool	76.7
			-110.0
Valco	52.1	Cryo trap heat	39.2
	50.0		35.0
GC 1	86.7	Cryo trap cool	33.9
	20.0		35.0
GC 2	25.4	Transfer line	78.6
	25.0		80.0
Trap compartment	35.1	MS transfer lines	130.1
	35.0		125.0
		Trap compartment airspace	37.3

Chromatogram - Flight 10 Boundary Layer





Problems

- Ozone interference from valve rotor

Compounds affected -	Acetaldehyde Propanal Butanal Pentanal	Acetone
----------------------	---	---------

Compounds NOT affected -	Methanol Ethanol	Methyl Ethyl Ketone 2-Pentanone
---------------------------------	---------------------	------------------------------------

- Pressure effects from undersized sampling pump (at >10 km)

Compounds affected -	Acetaldehyde Propanal Butanal Pentanal
----------------------	---

Compounds NOT affected -	Methanol Ethanol	Acetone Methyl Ethyl Ketone 2-Pentanone
---------------------------------	---------------------	---




U.S. PROTECTIVE G. ANALYST
Model 400






120V
1000W

GC-MS
Model 400

Model 400

NCAR Calibration Facility





Calibration Gas Intercomparison




HS Calibration Cylinder - ~ 1 ppm Acetone, Methanol, Ethanol, Acetaldehyde

FGCMS Calibration Cylinder - 556 ppb Acetone, 52 ppb MEK, 412 ppb Methanol, 54 ppb Ethanol, 102 ppb Acetaldehyde, 43 ppb Propanal, 11 ppb Butanal, 39 ppb Carbon Tetrachloride


Diluted into ultraclean zero air at ambient humidity

Compound	FGCMS standard (RF) (ppt/count)	HS standard (RF) (ppt/count)	% difference
Acetone	2.25×10^{-3}	2.0×10^{-3}	11.8
Methyl Ethyl Ketone	1.85×10^{-3}	n/a	
Methanol	2.35×10^{-2}	1.77×10^{-2}	28.1
Ethanol	1.97×10^{-2}	1.18×10^{-2}	50.3
Acetaldehyde	1.30×10^{-2}	1.07×10^{-2}	19.4
Propanal	3.40×10^{-3}	n/a	
Butanal	2.39×10^{-3}	n/a	


Sunrise Experiment



Flight 18 (Yakota - Kona Transit)




Flight 13: Yokota Local 1



Detection of Methyl Ethyl Ketone in the upper troposphere

NASA TRACE -P Experiment, Spring 2001, Missions 10 -20, Preliminary Data



Eric Apel, NCAR
Dan Riemer, UM

Potentially contributes to HO_x budget through photolysis
Lifetime approximately 1/4 that of acetone
Previously not included in budgets

Conclusions and Outlook:

Alcohols and carbonyls measured

- some problems w/aldehydes discovered post-flight
- in process of correcting data
- instrument being modified and improved

Unanswered questions: Ocean as a source of aldehydes?
comparison w/Singh after data correction

Potential Papers:

Instrument paper. The instrument used has never been described – it will be described in this paper

Distributions of methanol and ethanol during TRACE-P. A paper examining the distributions, sources, and sinks for these compounds.

Distributions of ketones during TRACE-P. A paper examining the distributions, sources, and sinks for these compounds.